



excelUNFORS Dosimeter Software


Create a folder intended to store the software:




Extract the files to the folder:



Open the install.vbs
file:




Install the downloaded autoit-v3-
setup.exe:







Click Finish:




Open excelUNFORS.au3




Choose the Dosimeter serial port. Please note clicking this dropdown scans all the cells in the current spreadsheet.





Click Browse for Template Folder to Choose where the Excel Templates are stored.




Click
Select



Click the Dropdown and Choose the Desired Template:



Type a value in the input box and press ENTER or click the Next Cell Arrow:



Click the kV button to input that value. Sets of values can be cycled through using the leftmost arrows.

